



## **GAZELLE** MELROSE TENNIS WHEELCHAIR

So you need to get to every ball, you need a stable platform to make shots from, and you need your chair to respond as if it was an extension of your body.

The new Gazelle has a light-weight frame designed to increase first push response. It has four contact points on court for a very stable ride, and because its custom built for you, it fits like a glove.

Built from aluminium or titanium your computer aided designed frame will offer maximum rigidity. This directs all kinetic pushing forces directly to the wheels with no energy lost in frame flexibility.

Custom built and with no adjustment matrix's, the Melrose Gazelle is a simple, uncluttered, very light and durable Tennis Wheelchair.

Get to the ball in a Melrose Gazelle.



3 wheel



4 wheel

## GAZELLE

### STANDARD SPECIFICATIONS

- Lightweight Aluminium 3 or 4 Wheel Frame
- Complete custom fit - built to individual specifications for Size, Rake and Camber
- Camber 10° - 20°
- Quick Release 24", 25" or 26" Spoked Wheels on Alloy Rims
- Smooth High Pressure Tyres
- Rollerblade Castors
- Black Fabric Sling and Back
- 5 Year Frame Warranty

### OPTIONS AND ACCESSORIES

- Super Lightweight Titanium 3 or 4 Wheel Frame
- Tapered Seat Width and Back
- Double or Single Swivel Anti-Tips
- Adjustable Back Upholstery
- Skirt Guards: Fabric Skirt Guards  
Alloy or Carbon Fibre Skirt Guards
- Straps and Bindings for Hips and Legs
- Cushion
- Smooth, High Pressure, Asymmetric Kenda Kaliente or Schwalbe Speedair Tyres
- Super light weight Spinerger Wheels - 24", 25" or 26" with colour option Carbon Cable Spokes
- Wheel Bag to accommodate four 24", 25" or 26" Wheels