



Do you need a chair that is built to take the hits and keeps coming back for more?

We use the highest quality materials and components to build our wheelchairs, many of the components are built in-house by our team of qualified engineers using the latest CNC technology. We also employ highly experienced fabricators who are working for a company that has over 24 years involvement in the sport.

As a small customer focused company we encourage you to talk directly with the fabricator of your new wheelchair to ensure you get the product you desire, this is a service we will continue to offer you.

If you are committed to the game, we are committed to providing you with the best equipment.

  
**Melrose**  
KIWI CONCEPT CHAIRS USA



Low Point Rhino (Defensive)



Mid Point Rhino



High Point Rhino (Offensive)

## RHINO

### STANDARD SPECIFICATIONS

- Fully custom built Offensive or Defensive Wheelchairs
- Choice of Pick Bar styles in Alloy
- Fixed Camber Bar between 12° and 22°
- Cross Spoked Wheels 24", 25" or 26" with High Pressure Tyres
- Hard Anodised Pushrims and Quick Release Axles
- Rollerblade Castors
- Adjustable Height Footplate
- Adjustable Height Back
- Adjustable Tension Seat and Back Upholstery, Fabric Skirt Guards, Foot Strap
- Single Wing (on Offensive only)

### OPTIONS AND ACCESSORIES

- Individual Designs for Wings can be accommodated - Double, Single and/or Inverted
- Tapered Seat Width and Back
- Ergo Seat and Back
- Double or Single Swivel Anti-Tips with Alloy or Plastic Plating Protection
- ABS or Alloy Footplate fixed or fully adjustable
- Offensive or Defensive style Steel Bolt on Pick Bar
- Spinergy Wheels SL or SLX with 1/2" or 5/8" Axles - 24", 25" or 26"
- Various Spoke Guard designs: -
  - Polycarb Offensive and Defensive
  - Polyethylene Offensive
  - Polyethylene / Red Rubber Defensive
  - Alloy Spoke Guards
- Red or Black Rubber coated Pushrims. Oversize Pushrims with 6 or 12 tabs
- Hip Strap, Back Bar Cover, Cushions 1" 2" 3" or 4"
- Velcro attached Chest Strap
- Alloy or ABS Skirt Guards
- Wheel Bag to accommodate four 24", 25" or 26" Wheels



MELROSE KIWI CONCEPT CHAIRS USA INC.

5671 E. Fountain Way • Fresno, CA 93727 • United States of America

Tel 866 764 5526 Tel (559) 348 2260 Fax (559) 292 0400 Email sales@melrosewheelchairs.com Web www.melrosewheelchairs.com