

RHINO

RUGBY SPORTS WHEELCHAIR



Do you need a chair that is built to take the hits and keeps coming back for more?

We use the highest quality materials and components to build our wheelchairs, many of the components are built in-house by our team of qualified engineers using the latest CNC technology. We also employ highly experienced fabricators who are working for a company that has over 26 years involvement in the sport.

As a small customer focused company we encourage you to talk directly with the fabricator of your new wheelchair to ensure you get the product you desire, this is a service we will continue to offer you.

If you are committed to the game, we are committed to providing you with the best equipment.



RHINO

RUGBY SPORTS WHEELCHAIR

STANDARD SPECIFICATIONS

ALUMINIUM FRAME

Fully custom built Offensive or Defensive Wheelchairs

Choice of Pick Bar styles in Alloy

Fixed Camber Bar between 12° and 22°

Cross Spoked Wheels 24", 25" or 26" with High Pressure Tyres

Hard Anodised Pushrims and Quick Release Axles

Rollerblade Castors

Adjustable Height Footplate

Adjustable Height Back

Adjustable Tension Seat and Back Upholstery, Fabric Skirt Guards, Foot Strap

Single Wing (on Offensive only)

OPTIONS AND ACCESSORIES

HEAT TREATED ALUMINIUM FRAME

TITANIUM FRAME

Individual Designs for Wings can be accommodated - Double, Single and/or Inverted

Choice of Pick Bar styles in Steel or Titanium

Tapered Seat Width and Back

Ergo Seat and Back

Double or Single Swivel Anti-Tips with Alloy or Plastic Plating Protection

ABS or Alloy Footplate fixed or fully adjustable

Offensive or Defensive style Steel Bolt on Pick Bar

Spinergy Wheels SL or SLX with 1/2" or 5/8" Axles - 24", 25" or 26"

Various Spoke Guard designs: -

- Polycarb Offensive and Defensive
- Polyethylene Offensive
- Polyethylene / Red Rubber Defensive
- Alloy Spoke Guards

Red or Black Rubber coated Pushrims. Oversize Pushrims with 6 or 12 tabs

Hip Strap, Back Bar Cover, Cushions 1" 2" 3" or 4"

Velcro attached Chest Strap

Alloy or ABS Skirt Guards

Wheel Bag to accommodate four 24", 25" or 26" Wheels

