

RUGBY LEAGUE SPORTS WHEELCHAIR

# WARRIOR



**MelroseHeat**   
HEAT TREATED ALUMINIUM FRAME



As a recent sport, the technology for Rugby League wheelchairs has already been developed with the light weight of the Melrose Basketball wheelchair and the strength of the Melrose Rugby wheelchair. We have combined these two ideas with a light weight 1" frame built with 1.5" bottom rails. The bottom rails are further reinforced with 1.5" tubes to the capped 1.5" wing for additional strength. Although Wheelchair Rugby League is a non

contact sport, experience has taught us that a strong wheelchair will perform better over extended time.

All our Rugby wheelchairs are built by highly experienced fabricators who understand the sport of rugby. We use the highest quality materials and components to build our wheelchairs.



# WARRIOR

## RUGBY LEAGUE SPORTS WHEELCHAIR

### STANDARD SPECIFICATIONS

Complete custom fit - built to individual specifications for size, rake and camber

Lightweight Aluminium frame

Fixed camber bar between 12° and 20°

24", 25" or 26" spoked wheels with high pressure tyres on alloy rims

Hard anodised pushrims and quick release axles

Rollerblade castors

Adjustable height footplate

Cushion

Adjustable height back

Adjustable tension seat sling, back upholstery, fabric skirt guards, foot strap

### OPTIONS AND ACCESSORIES

#### HEAT TREATED ALUMINIUM FRAME

Super lightweight Titanium frame

Tapered seat width and back

Adjustable tension back upholstery

Ergonomically designed seating system (Ergo Seat)

Hip strap, back bar cover, cushions 1" 2" 3" or 4"

Double or single swivel anti-tips

Aluminium, Alloy or ABS skirt guards

Rubber coated Hand Rims - Red or Black

Spoke guards - clear or black

Tyres, smooth, high pressure IRC and Kenda Iron Cap

Super light weight Spinergy wheels - 24", 25", 26" or 700c with colour option carbon cable spokes

ABS or alloy footplate fixed or fully adjustable

Various spoke guard designs:

- Polycarb

- Polyethylene

- Polyethylene / red rubber

- Alloy spoke guards

Wheel bag to accommodate (two or four) 24", 25" or 26" wheels

